**Instruction Sheet for testing of any combination of the following drugs:
OPI/MDMA/COC/AMP/THC**

A rapid test for the simultaneous, qualitative detection of multiple drugs in human urine.
For self-testing use.

【INTENDED USE AND SUMMARY】

The Multi 5 Drugs Rapid Test Panel is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in urine at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
Opiate (OPI)	Morphine	2000
Methylenedioxy-methamphetamine (MDMA)	d,l-Methylenedioxy-methamphetamine	500
Cocaine (COC)	Benzoylcegonine	300
Amphetamine (AMP)	d-Amphetamine	500
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	50

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

【PRINCIPLE】

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test region of the specific drug dipstick. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test region.

A drug-positive urine specimen will not generate a colored line in the specific test region of the dipstick because of drug competition, while a drug-negative urine specimen will generate a line in the test region because of the absence of drug competition.

To serve as a procedural control, a colored line will always appear at the control region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

【PRECAUTIONS】

- For self-testing use.
- The Test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

【STORAGE AND STABILITY】

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The Test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【SPECIMEN COLLECTION AND PREPARATION】**Urine Assay**

The urine specimen should be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be filtered or allowed to settle to obtain a clear specimen for testing. Avoid consuming alcohol in connection with the test, as this may give a positive result.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing.

【MATERIALS】**Materials Provided**

- Test Panel

- Package Insert

Materials Required But Not Provided

- Timer

【DIRECTIONS FOR USE】


Allow the test and urine specimen to reach room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test panel from the sealed pouch and use it within one hour.
- Remove the cap.
- With the arrow pointing toward the urine specimen, immerse the test panel vertically in the urine specimen for at least 10 to 15 seconds. **Immerse the dipstick to at least the level of the wavy lines, but not above the arrow on the test panel.**
- Replace the cap and place the test panel on a non-absorbent flat surface.
- Start the timer and wait for the colored line(s) to appear.
- The drug strip result should be read at **5 minutes**. Do not interpret the result after 10


minutes.

TEST PROCEDURE Panel

- Remove the cap.
- Do not immerse past the max line
- Immerse 10-15 seconds

Read Drug Result at 5 minutes.

**【INTERPRETATION OF RESULTS】**

(Please refer to the illustration above)

NEGATIVE: A colored line appears in the control region (C) and colored line appears in the test region (T). This negative result means that the concentrations in the urine sample are below the designated cut-off levels for a particular drug tested.

***NOTE:** The shade of the colored lines(s) in the test region (T) may vary. The result should be considered negative whenever there is even a faint line.

POSITIVE: A colored line appears in the control region (C) and no line appears in the test region (T). The positive result means that the drug concentration in the urine sample is greater than the designated cut-off for a specific drug.

INVALID: No line appears in the control region (C). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Read the directions again and repeat the test with a new test. If the result is still invalid, contact your manufacturer.

【QUALITY CONTROL】

A procedural control is included in the test. A line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

- The Multi 5 Drugs Rapid Test Panel provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography /mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may cause erroneous results.
- A positive result does not indicate level or intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- This test does not distinguish between drugs of abuse and certain medications.
- A positive test result may be obtained from certain foods or food supplements.

【PERFORMANCE CHARACTERISTICS】**Accuracy
% Agreement with GC/MS**

	OPI 2,000	MDMA 500	COC 300	AMP 500	THC 50
Positive Agreement	96.7%	98.1%	98.2%	99.1%	97.9%
Negative Agreement	93.8%	99.3%	97.8%	98.6%	98.1%
Total Results	95.2%	98.8%	98.0%	98.8%	98.0%

% Agreement with Commercial Kit

	OPI 2,000	MDMA 500	COC 300	AMP 500	THC 50
Positive Agreement	*	>99.9%	>99.9%	>99.9%	>99.9%
Negative Agreement	*	>99.9%	>99.9%	>99.9%	>99.9%
Total Results	*	>99.9%	>99.9%	>99.9%	>99.9%

*Note: Based on GC/MS data instead of Commercial Kit.

Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical card of coded specimens, containing drugs at concentrations of negative, $\pm 50\%$ and $\pm 25\%$ cut-off level, was labeled, blinded and tested at each site. The results gained $\geq 75\%$ accuracy in $\pm 25\%$ cut-off level specimen and 100% accuracy in negative and

$\pm 50\%$ cut-off level specimen.

Analytical Sensitivity

A drug-free urine pool was spiked with drugs at the listed concentrations. The results are summarized below.

Drug Concentration Cut-off	OPI 2,000		MDMA 500		COC 300		AMP 500		THC 50	
Range	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	27	3	25	5	26	4	25	5	26	4
Cut-off	15	15	14	16	13	17	15	15	14	16
+25% Cut-off	5	25	4	26	3	27	3	27	3	27
+50% Cut-off	0	30	0	30	0	30	0	30	0	30
+300% Cut-off	0	30	0	30	0	30	0	30	0	30

Analytical Specificity

The following table lists the concentrations of compounds (ng/mL) that are detected as positive in urine by the Multi 5 Drugs Rapid Test Panel at 5 minutes.

Analytes	conc. (ng/mL)	Analytes	conc. (ng/mL)
MORPHINE/OPIATE (OPI 2,000)			
Codeine	2,000	Morphine	2,000
Ethylmorphine	3,000	Norcodeine	25,000
Hydrocodone	50,000	Normorphine	50,000
Hydromorphone	15,000	Oxycodone	25,000
Levorphanol	25,000	Oxymorphone	25,000
6-Monoacetylmorphine	3,000	Procaine	50,000
Morphine 3- β -D-glucuronide	2,000	Thebaine	25,000
METHYLENEDIOXYMETHAMPHETAMINE (MDMA 500) Ecstasy			
($\pm$) 3,4-Methylenedioxy-methamphetamine HCl	500	3,4-Methylenedioxyethyl-amphetamine	300
($\pm$) 3,4-Methylenedioxyamphetamine HCl	3,000		
COCAINE (COC 300)			
Benzoylcegonine	300	Cocacethylene	20,000
Cocaine HCl	200	Ecgonine	30,000
AMPHETAMINE (AMP 500)			
D,L-Amphetamine sulfate	150	Phentermine	500
L-Amphetamine	12,500	Maprotiline	25,000
($\pm$) 3,4-Methylenedioxy amphetamine	250	Methoxyphenamine	3,000
		D-Amphetamine	500
MARIJUANA (THC 50)			
Cannabinol	35,000	Δ^8 -THC	17,000
11-nor- Δ^8 -THC-9 COOH	30	Δ^9 -THC	17,000
11-nor- Δ^9 -THC-9 COOH	50		

Fifteen (15) urine samples of normal, high, and low specific gravity ranges (1.005-1.045) were spiked with drugs at 50% below and 50% above cut-off levels respectively. The Multi 5 Drugs Rapid Test Panel was tested in duplicate using fifteen drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with drugs at 50% below and 50% above cut-off levels. The spiked, pH-adjusted urine was tested with the Multi 5 Drugs Rapid Test Panel. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug positive urine containing above related calibrator substances. The following compounds show no cross-reactivity when tested with the Multi 5 Drugs Rapid Test Panel at a concentration of 100 μ g/mL.

Non Cross-Reacting Compounds

Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labetalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Isoxsuprine	Sulfamethazine
l-Ascorbic acid	Digoxin	d,l-Propanolol	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,

Atropineβ-EstradiolNiacinamide3-acetateBenzilic acidEstrone-3-sulfateNifedipineTetrahydrocortisoneBenzoic acidErythromycinNorethindroneTetrahydrozolineBilirubinFenoprofenNoscapineThiamined,l-BrompheniramineFurosemide d,l-OctopamineThioridazineGentisic acidOxalic acid d,l-TyrosineHemoglobinOxolinic acidTolbutamideChloral hydrateHydralazineOxymetazolineTriamtereneChloramphenicolHydrochlorothiazidePapaverineTrifluoperazineChlorothiazideHydrocortisonePenicillin-GTrimethoprimd,l-Chlorpheniramineo-Hydroxyhippuric acidPerphenazine d,l-TryptophanChlorpromazine3-HydroxytyraminePhenelzineUric acidCholesterol d,l-IsoproterenolPrednisoneVerapamilClonidine

【BIBLIOGRAPHY】

1. Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

2. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 6th Ed. Biomedical Publ., Foster City, CA 2002.

Index of Symbols			
	Consult instructions for use		Contains sufficient for <n> test
	Do not use if package is damaged and consult instructions for use		Catalogue number
	Manufacturer		Use-by date
	Store between 2-30 °C		Batch code
	Do not reuse		Authorized Representative

Hangzhou AllTest Biotech Co.,Ltd.
#550,Yinhuai Street
Hangzhou Economic & Technological Development Area
Hangzhou, 310018 P.R. China
Web: www.alltests.com.cn Email: info@alltests.com.cn

MedUnion S.L.
Carrer de Tapioles, 33, 2-1, 08004,
Barcelona, Spain
Email: rep@themedunion.com

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